

Product Information Sheet



Product Code: 50-001-0185-0637

Product Name: NEANN Patient Transport Kit

Purpose / Function / Design

Easy to use, well-padded bag, designed for Intensive Care Units or Hospital in the Home to store and transport Medical equipment and accessories

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Clear windows on pockets and side of bag: Allows ease of adding identification or Logo's
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling: Allows ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Magnetised, replaceable and padded aVector handles: Minimizes cross contamination
- Padded, adjustable and replaceable shoulder sling
- Fidlock buckles on sling: quick release in emergency situation and ease of reattachment
- Magnetic and replaceable sharps holder with locator : bag supplied with 600ml holder, option to purchase Larger Bio-can size
- Foam rigidity to outer walls, lid and base : maintains form and adds impact resistance
- Heavy duty zip placed to give maximum access to interior
- Single zip pull handle : ease of use if hands are gloved, one hand operation
- "Sharps attach" print : quickly re-attach Sharps holder, no search time
- Three large external pockets with Magnetic pocket flaps: quick access and self closing

Internal

- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Side pocket windows: Improves on scene efficiencies through labelling of contents, personalized by end user
- Five stash pockets and one sleeve in main compartment
- Two pockets inside lid
- Main compartment dimensions will fit $\frac{3}{4}$ Airways and $\frac{3}{4}$ Drug Kits (ordered separately)

Specifications

External dimensions:	450mm (L) x 260mm (W) x 235mm (H)
Total volume:	26 Litres
Weight when empty:	1.6kg
Maximum recommended carry weight:	7kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated fabric
Internal structure:	Foam

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Information Sheet



Optional purchase is 1.4L Bio-can holder 50-001-8082-0386 , as pictured

