

Product Information Sheet



Product Code: 50-001-0630-0192

Product Name: NEANN PRK

Purpose / Function / Design

Designed to store and transport Primary Response equipment, contents can be Customised to be either Airway / IV / or Trauma focused

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Kit clearly named: Allows ease of identification
- Clear ID window: customize bag with business card
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling: Allows ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Reinforced base: Abrasion resistant fabric
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Magnetised, replaceable and padded aVector handles: Minimizes cross contamination
- Padded, adjustable and replaceable shoulder sling
- Padded, adjustable backpack straps with cover : keep hands free when carrying and protects straps when stowed
- Side release buckles on backpack straps: quick release in emergency situation and ease of reattachment
- Magnetic and replaceable sharps holder
- “Sharps attach” print on bag: quickly re-attach Sharps holder, no search time
- Foam rigidity to outer walls, lid and base : maintains form and adds impact resistance
- Heavy duty zip placed to give maximum access to interior
- Yellow zip pull tab extension: good visibility and ease of grab when wearing gloves
- Lockable zip on main compartment : extra security of contents
- Large top pocket : holds drug or airways kits (separate products)
- Cervical collar holder
- Smaller side pocket for quick access items

Internal

- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Light coloured fabric: Allows contamination to be easily identified
- Magnetic internal lid surface : can customize bag by adding either Trauma, Airways or IV panel (separate products : NEANN PRK (Trauma panel) 50-001-0634-0422, NEANN PRK (Airway panel) 50-001-0636-0422 , NEANN PRK (IV panel) 50-001-0638-0422) ease of cleaning
- Clear sleeve in lid pocket: Provides easy storage of documents
- Lid pocket can fit NEANN Advanced Drug Kit 50-001-0635-0386 (separate product)
- Dividers in main compartment : organize contents, dividers are removable

Product Information Sheet



Specifications

External dimensions:	570mm (L) x 300mm (W) x 300mm (H)
Total volume:	51 Litres
Weight when empty:	2.7kg
External fabric:	Medium weight reinforced EnvironFrei
Staerksager:	Abrasion resistant fabric
Internal fabric:	Lightweight coated fabric
Internal structure:	Recycled Corflute, Foam



Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365