

Product Information Sheet



Product Code: 50-001-8728-0369

Product Name: NEANN EZ-IO Pouch

Purpose / Function / Design

Designed to store and transport an intraosseous infusion device and five needle sets.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing grab handle : Allows ease of cleaning and minimizes cross contamination
- Kit clearly named: Allows ease of identification
- Foam rigidity to top and bottom for maintaining form and impact resistance
- Lockable zip : security of contents

Internal

- MEDIS labelling: Improves on scene efficiencies through quick identification of contents
- Customized design: Eliminates “Search & Select” issues
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified

Specifications

External dimensions:	140mm (L) x 120mm (W) x 100mm (H)
Total volume:	1.6 Litres
Weight when empty:	182 g
Maximum recommended carry weight:	590 g
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Fabric
Internal structure:	Foam

Date: March 2022

Product Information Sheet



Pouch supplied empty

Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365