



Buying Guide

Compare Veinlite® Models by Application

	Portable LED Devices					
	EMS PRO	LED+	LEDX	PEDI2	NEO	R
Adult IV Access	✓	✓				
Pediatric IV Access	✓	✓		✓		
Emergency	✓	✓				
Neonatology					✓	
Sclerotherapy & Phlebectomy			✓			
Oncology & Infusion		✓				
Radiology & Nuclear Medicine		✓				
Home Health		✓				
Research						✓

Compare Portable Devices by Features

	Portable LED Devices					
	EMS PRO	LED+	LEDX	PEDI2	NEO	R
Number of LEDs	24	28	32	14	3	3
Color of LEDs	Orange (12) + Red (4) + White (8)	Orange (22) + Red (6)	Orange (24) + Red (8)	Orange (10) + Red (4)	Orange (veins) + Green (arteries) + White (exam)	Orange (veins) + Green (arteries) + White (exam)
Color Selection	Yes	Yes	Yes	Yes	Yes	Yes
Brightness Selection	No	No	No	Yes	Yes, 3 Levels	Yes, 3 Levels
Ring Diameter (mm)	21	27	31	18	N/A	N/A
Battery	2 x AA Alkaline	Rechargeable Lithium Ion	Rechargeable Lithium Ion	2x AAA Alkaline	1 x AAA Alkaline	1 x AAA Alkaline
Battery Run Time (mins)	550	350	130 / 160 (32 / 24)	-	-	-
Dimensions (mm)	115 x 55 x 21	99 x 57 x 21	102 x 65 x 21	85 x 40 x 17	96.1 x 19 x 16.5	96.1 x 19 x 16.5
Weight, inc. battery (g)	108	77	83	64	29	29