

# Product Information Sheet



**Product Code:** 50-001-0627-0386

**Product Name:** Neann 3/4 Drug Kit

## **Purpose / Function / Design**

Designed to store and transport IV / Cannulation gear, drugs ampoules, bottles and blister packs

## **Features & Benefits**

### **External**

- Reflective tape: Provides high visibility in low light conditions
- Kit clearly named: Allows ease of identification
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- Heat transfer artwork: Allows ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing handle: Allows ease of cleaning and minimizes cross contamination
- Rigidity to lid and base : maintains form and adds impact resistance
- Heavy duty lockable zip : maximum access to interior and security
- Yellow zip pull tab extension: good visibility and ease of grab when wearing gloves
- Clear sleeve: Provides easy storage of documents

### **Internal**

- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Internal flip panel : increased storage area
- Internal flip panel has colour coded webbing attach: keeps layout consistent across fleet
- Elasticated pockets with securing strap: larger items/ blister packs
- Removable Drug and IV Cannula loops : ease of re-stocking and content Layout can be customized
- Potential to order Medis printed drug loops (sold separately)

## **Specifications**

External dimensions:	350mm (L) x 215mm (W) x 100mm (H)
Total volume:	7.5 Litres
Weight empty:	867g
Maximum recommended carry weight:	3kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Polyester
Internal structure:	Recycled Corflute



# Product Information Sheet



Pan a and b - Cannulation



Pan c and d – Drugs



Head Office  
RAPP Australia Pty. Ltd.  
160 McClelland Ave, PO Box 120,  
Lara, Victoria, 3212, Australia  
ABN 35 054 189 365

# Product Information Sheet



- 2 x Tourniquet
- 1 x Shaver
- 10 x Swabs
- 12 x IV Cannula
- 5 x IV Dressing
- 4 x Reflux valves/bungs
- 2 x 3 way tap
- 1 x Tape
  
- 1 x IN Syringe
- 1 x ampoule breaker
  
- 5 x Adhesive dressing Large/IO needle
- Needles (up to x18)
- Syringes (1ml x3, 3ml x3, 5ml x6, 10ml x4, 20ml x1)

# Product Information Sheet



Contents (from Left to right)

- Sodium Chloride polyamps x6 ,  
Metoclopramide x2
- Atropine x4, Salbutamol x6
- Blister pack holder x4 (Paracetamol,  
Ibuprofen, Aspirin, Fexofenadine)
- 
- adrenaline x8, naloxone x6
- midazolam x2, midaz 3ml x2,  
amiodarone x3, Adrenaline 1:10,000 x2,  
Lignocaine 2% x2
- adenosine x5, benzylpenicillin x2,  
hydrocortisone x1
- GTN spray, Epipen, glucose paste x2,  
glucagon x1, methoxyflurane box

