

Allied LSP - Brass Laser Spectrum Oxygen Regulators

Model Numbers and Options

Part Number	0-25 Regulator 0, 0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20 and 25 Litres per Minute	0-8-15 Regulator 0, 8 and 15 Litres per Minute	Australian Standard Diameter Indexed Low Pressure Outlets (400kPa)
50-002-32292550ASTL25-0000	Yes		2
50-002-32292502AU-0000		Yes	0
50-002-32292502AU2-0000		Yes	2

If you cannot see a model that will meet your needs,
please contact us and we will probably be able to build it up for you.