

Product Information Sheet



Product Code: 50-001-0128-0192

Product Name: Hospital Response Kit

Purpose / Function / Design

Designed for hospitals to store and transport emergency response equipment

Features & Benefits

External

- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- aVector logo : ease of cleaning and minimizes cross contamination
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Magnetised, replaceable and padded aVector handles: Minimizes cross contamination
- Padded, adjustable and replaceable shoulder sling
- Foam rigidity to outer walls, lid and base : maintains form and adds impact resistance
- Heavy duty zip: placed to give maximum access to interior
- Clear ID window: bag customization option with business card

Internal

- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Clear pocket: Provides visibility of contents
- Magnetic flap closure on pouch : one hand usage and auto close
- Divider separated compartment: contents can drop in but still be visible
- Dividers are movable, removable and wipeable: Allows ease of cleaning and minimizes cross contamination

Specifications

External dimensions:	380mm (L) x 260mm (W) x 150mm (H)
Total volume:	14.5 Litres
Weight when empty:	1.25kg
Maximum recommended carry weight:	6kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute & Foam

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365