

Product Information Sheet



Product Code: 50-001-0634-0422

Product Name: NEANN PRK (Trauma Panel)

Purpose / Function / Design


Designed for the NEANN Primary Response Kit to store and transport trauma equipment.

Features & Benefits

- Magnetic removable panel: ease of cleaning and allows the panels to be replaced or updated.
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- MEDIS labelling: Improves on scene efficiencies through quick identification of contents
- Clear pouches: easily view and secure small contents
- Magnetic lid flap on pouches : one hand usage, auto close

Specifications

Dimensions:	500mm (L) x 270mm (W)
Weight when empty:	0.39kg
Stocked weight:	1.13kg
Fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute



Product Information Sheet



Shown in NEANN PRK 50-001-0630-0192
(Separate Product)

Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365