

Product Information Sheet



Product Code: 50-001-0636-0422

Product Name: NEANN PRK (Airway Panel)

Purpose / Function / Design

Designed for the NEANN Primary Response Kit to store and transport airway equipment.

Features & Benefits

- Magnetic removable panel: ease of cleaning and allows the panels to be replaced or updated.
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- MEDIS labelling: Improves on scene efficiencies through quick identification of contents

Specifications

Dimensions:	500mm (L) x 270mm (W)
Weight when empty:	0.36kg
Stocked weight:	1.23kg
Fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute



Product Information Sheet



Shown in NEANN PRK 50-001-0630-0192
(Separate Product)