

# Product Information Sheet



**Product Code:** 50-001-0639-0422

**Product Name:** NEANN TEK+ (Lid Insert Trauma)

## **Purpose / Function / Design**

Designed for the NEANN TEK + to store and transport trauma equipment.

## **Features & Benefits**

- Magnetic removable panel: ease of cleaning and allows the panels to be replaced or updated.
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking. Including a section of loops left unlabelled so it can be personalized by end user.
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination.
- Lightly coloured fabric: Allows contamination to be easily identified.
- MEDIS labelling: Improves on scene efficiencies through quick identification of contents.


## **Specifications**

Dimensions:	380mm (W) x 27mm (H)
Weight when empty:	228g
Stocked weight:	644g
External fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute

**Purchased by:**

**Date:**

## **Images:**



Shown in NEANN TEK+ 50-001-0631-0192  
(Separate Product)