



**Laerdal**  
helping save lives



# The Bag II

## Transparent Disposable Resuscitators

All of the outstanding features available in a reusable resuscitator with the convenience of a disposable resuscitator. Performance, versatility and convenience when you need it most.

In keeping with our commitment to provide quality products to help you save lives, Laerdal's The Bag II is an important supplement to our complete product offering for Prehospital, Hospital, and the Military.



# The BAG II Product Features and Benefits

## Three Bag Volumes; Adult, Child, Infant

- Appropriate volume for the patient
- Allows the clinician a closer match to the volume needs of the patient without hyperinflation

## Disposable

- Single use, helps reduce risk of cross contamination
- No loss of reusable inventory

## Latex-Free

- Safe for latex-sensitive patients and staff

## 35cm Pop Off Available

- Infant/Child models only

## Masks

- Soft cushion with clear housing to view patient's skin color

## Temperature Range

- Performs well in both high and low temperature conditions
- Operating temperature 14-122°F (-10 to +50°C)

## PVC Construction

- Superior feel of lung compliance
- Transparent construction
- Ergonomic design

## Adjustable Hook and Loop Handle

- Reduces hand fatigue

## Reinforced Non-Kink O<sub>2</sub> Tubing

- Prevents flow obstruction

## Oxygen Reservoir

- Allows 99-100% O<sub>2</sub> administration

## Universal Star-Lumen Oxygen Tubing

- No need to have adapter, tubing fits directly on flowmeter

## ORDERING INFORMATION:

### Singles

845131	The BAG II Resuscitator Infant w/mask #1
845123	The BAG II Resuscitator Child w/mask #2
845121	The BAG II Resuscitator Child w/mask #3
845141	The BAG II Resuscitator Adult w/mask #4
845111	The BAG II Resuscitator Adult w/mask #5
845151	The BAG II Resuscitator Adult w/mask #5 with inflation port
845152	The BAG II Resuscitator Child w/mask #3 with inflation port

### Packs of 12

845211	The BAG II Resuscitator Adult w/mask #5
845221	The BAG II Resuscitator Child w/mask #3
845223	The BAG II Resuscitator Child w/mask #2
845231	The BAG II Resuscitator Infant w/mask #1
845241	The BAG II Resuscitator Adult w/mask #4
845252	The BAG II Resuscitator Child mask #3 with inflation port

### Accessories

845240	Laerdal Disposable PEEP Valve Assembly, pkg 10
845255	Laerdal Disp. Mask #1 w/inflation port, pkg 20
845256	Laerdal Disp. Mask #2 w/inflation port, pkg 20
845250	Laerdal Disp. Mask #3 w/inflation port, pkg 20
845260	Laerdal Disp. Mask #4 w/inflation port, pkg 20
845270	Laerdal Disp. Mask #5 w/inflation port, pkg 20
850900	LSR Manometer connector

Note: All child and infant resuscitators include manometer port and 35 cm H<sub>2</sub>O pop-off.

## REGULATORY STANDARDS AND APPROVALS

Europe: CE marked as class IIa device in accordance with Council Directive 93/42/EEC Medical Device Directive, as amended by Council Directive 2007/47/EC.  
USA: Cleared under 510(k) as a class 2 device.  
Canada: Holds a class 2 medical device license



The product is in compliance with the essential requirements of Council Directive 93/42/EEC as amended by Council Directive 2007/47/EC.

Laerdal Medical is certified by DNV GL Presafe AS to ISO 13485:2016



**Laerdal**  
helping save lives