

Product Information Sheet



Product Code: 50-001-0170-0192

Product Name: HITH Vehicle Kit

Purpose / Function / Design

Designed with maximum Infection control finishes and materials for Hospital In The Home, to store various equipment. Contains cooler pouch.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling: Allows ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Heavy duty double zip with pull handle: placed to give maximum access to interior
- Padded top handles : grab point to remove bag from vehicle
- Padded, adjustable and replaceable shoulder sling
- Large side release buckle: secure trolley attachment
- Seat belt loops: safely secures load in vehicle
- External compartment with magnetic lid flap: auto close and security for contents

Internal

- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Removable and movable dividers
- End stash pocket
- Cooler pouch with clear lid and magnetic flap opening

Specifications

External dimensions:	380mm (L) x 320mm (W) x 340mm (H)
Total volume:	41 Litres
Weight when empty:	2.5 kg
Maximum recommended carry weight:	8 kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute and Foam

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Information Sheet



Additional Colourways:



50-001-0170-0369



50-001-0170-0637