

Product Information Sheet



Product Code: 50-001-0202-0192

Product Name: HITH Backpack

Purpose / Function / Design

Designed to store and transport hospital in the home equipment.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling: ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Magnetised, replaceable and padded aVector handles: Minimizes cross contamination
- Attachment point for padded, adjustable and replaceable shoulder sling – purchased separately if required (Padded Sling Straight Grey Golf Hook 98-001-8085-0422)
- Padded, adjustable backpack with sternum strap feature: keep hands free when carrying
- Attachment point for magnetic and replaceable sharps holder - purchased separately if required (code subject to size required)
- Foam rigidity to outer walls, lid and base: maintains form and adds impact resistance
- Heavy duty zip placed to give maximum access to interior
- Front external padded pocket - Magnetised flap closure

Internal

- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Removable panel work surface : clean prep surface for chosen location
- Internal lid has clear pockets with magnetic flaps: one hand usage and auto close
- Can fit pouches in base – purchased separately if required (Small Pouch - 50-001-0194-0422, Large Pouch - 50-001-0192-0422, HITH Observation Pouch (Red) - 50-001-0189-0192)

Specifications

External dimensions:	450mm (L) x 330mm (W) x 230mm (H)
Total volume:	34.2 Litres
Weight when empty:	1.88kg
Maximum recommended carry weight:	3kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Polyester
Internal structure:	Polyfoam & Recycled Corflute

Date: February 2023

Product Information Sheet



Photos:



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Information Sheet



With Padded Sling (Straight Grey Golf Hook) 98-001-8085-0422 (separate product):



With 600ml Sharps container holder (2 mags) 50-001-8088-0386 (separate product):



With pouches (separate products):



Small Pouch - 50-001-0194-0422

Large Pouch - 50-001-0192-0422

HITH Observation Pouch (Red) - 50-001-0189-0192