



Product Information Sheet

Product Code: 50-001-0309-0192

Product Name: NEANN Surf PRK

Purpose / Function / Design

Designed for Surf Life Saving as a Primary Response Kit, can store and transport various trauma and airways equipment

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Reinforced abrasion resistant base: increase the life of the bag
- Wipeable external fabric : Easy to clean and minimizes cross contamination
- aVector labelling: Allows ease of cleaning and minimizes cross contamination
- Padded carry handle for user comfort
- Padded, adjustable and replaceable shoulder sling
- Kit clearly named: Allows ease of identification
- External label windows : Provides flexibility with naming customization
- Heavy duty coated zips : resist corrosion and moisture
- Built in, abrasion resistant backpack strap cover Conceals backpack when not required reducing catch hazards
- Front zip pocket with multiple pockets and holders
- External side stash pocket

Internal

- Customized design : Eliminates “Search & Select” issues
- Everything has a place and is in its place : Eliminates overstocking providing a cost saving
- Water resistant fabric : Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric : High visibility, allows for contamination to be easily identified
- Clear windows on pockets : Provides flexibility with labelling
- Moveable dividers in Main compartment
- Elastic loops in varying sizes for varied usage

Specifications

External dimensions:	580mm (L) x 300mm (W) x 300mm (H)
Total volume:	52 Litres
Weight when empty:	3.13kg
Maximum recommended carry weight:	10kg
External fabric:	Medium weight reinforced VironFrei
Staerksager:	Abrasion resistant fabric
Internal fabric:	Lightweight coated fabric
Internal structure:	Recycled Corflute

Date: Jan 2020

