



Product Information Sheet

Product Code: 50-001-0329-0386

Product Name: Surf OTK

Purpose / Function / Design

Designed for Surf Life Saving NSW to store and transport Oxygen cylinder and related Oxygen Therapy equipment

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Reinforced abrasion resistant base: increase the life of the bag
- Wipeable external fabric : Easy to clean and minimizes cross contamination
- Heat transfer artwork: Allows ease of cleaning and minimizes cross contamination
- Padded carry handle for user comfort
- Padded, adjustable and replaceable shoulder sling
- Kit clearly named: Allows ease of identification
- External label windows : Provides flexibility with naming
- Heavy duty coated zips : resist corrosion and moisture
- Built in, abrasion resistant backpack strap cover Conceals backpack when not required reducing catch hazards
- Front drop in pocket with side release clip : Provides flexibility in allowing for additional products to be stored. (holds Surf Defib Case 50-001-0304-0192)
- External access to Oxygen regulator

Internal

- Customized design : Eliminates “Search & Select” issues
- Everything has a place and is in its place : Eliminates overstocking providing a cost saving
- Water resistant fabric : Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric : High visibility, allows for contamination to be easily identified
- Clear windows on pockets : Provides flexibility with labelling
- Clear sleeve : Allows for layout guide
- End pocket to hold suction and canister : Contains potential spill, reducing cross contamination

Specifications

External dimensions:	600mm (L) x 300mm (W) x 300mm (H)
Total volume:	54 Litres
Weight when empty:	3.19kg
Maximum recommended carry weight:	10kg
External fabric:	Medium weight reinforced VironFrei
Staerksager:	Abrasion resistant fabric
Internal fabric:	Lightweight coated Polyester
Internal structure:	Recycled Corflute

Purchased by: Surf NSW

Date: 22.3.19



