

# Product Information Sheet



**Product Code:** 50-001-0649-0192

**Product Name:** NEANN Transport Kit+

## **Purpose / Function / Design**

Easy to use, well-padded bag, designed for Intensive Care Units or Hospital in the Home to store and transport Medical equipment and accessories

## **Features & Benefits**

### **External**

- Reflective tape: Provides high visibility in low light conditions
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling: Allows ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Magnetised, replaceable and padded aVector handles: Minimizes cross contamination
- Foam rigidity to outer walls, lid and base : maintains form and adds impact resistance
- Heavy duty double zip with pull handle: placed to give maximum access to interior
- Long pull tabs for ease of use if hands are gloved
- Eight large external pockets with Magnetic pocket flaps: quick access and self closing (each 14cm W x 18cm H x 8cm D )
- Clear windows on pockets: Allows ease of adding identification or Logos
- Grab Handle on front: ease of moving the bag
- Security loop: lock main compartment with security tag
- Built in: "Sharps attach" locator for magnetic sharps holder: quickly re-attach Sharps holder, no search time. Option to purchase 600ml or 1.4L Bio-Can Sharps Container Holder (Separate products, not included in bag)

### **Internal**

- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Side pocket windows: Improves on scene efficiencies through labelling of contents, personalized by end user
- Five stash pockets and one sleeve in main compartment
- Two pockets inside lid
- Main compartment dimensions 63cm (L) x 27cm (W) x 27.5cm (H) will fit Neann Advanced Drug Kit - 50-001-0635-0386 / Neann  $\frac{3}{4}$  Airways Kit - 50-001-0633-0422 / Neann  $\frac{3}{4}$  Drug Kit 50-001-0627-0386 (ordered seperately).

# Product Information Sheet



## Specifications

|                                   |                                      |
|-----------------------------------|--------------------------------------|
| External dimensions:              | 630mm (L) x 430mm (W) x 275mm (H)    |
| Total volume:                     | 74 Litres                            |
| Weight when empty:                | 2.73kg                               |
| Maximum recommended carry weight: | 8kg                                  |
| External fabric:                  | Medium weight reinforced EnvironFrei |
| Internal fabric:                  | Lightweight coated Polyester         |
| Internal structure:               | 12mm Polyfoam                        |



# Product Information Sheet



Optional purchase is 600ml Sharps Container Holder 50-001-3079-0386 (as pictured) or Neann Bio-can 1.4L Sharps Container Holder 50-001-8184-0386

(separate products, not included in the bag)



Head Office  
RAPP Australia Pty. Ltd.  
160 McClelland Ave, PO Box 120,  
Lara, Victoria, 3212, Australia  
ABN 35 054 189 365