

Product Information Sheet



Product Code: 50-001-0633-0369

Product Name: NEANN ¾ Airway Kit

Purpose / Function / Design

Designed to store and transport airways access equipment.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Kit clearly named: Allows ease of identification
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling : ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing Handle: Allows ease of cleaning and minimizes cross contamination
- Yellow zip pull tab extension: good visibility and ease of grab when wearing gloves
- Lockable zip : extra security of contents
- Clear sleeve: Provides easy storage of documents

Internal

- MEDIS labelling: Improves on scene efficiencies through quick identification of contents
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Connection point for Paediatric airway panel (separate product 50-001-0642-0422)

Specifications

External dimensions:	350mm (L) x 210mm (W) x 100mm (H)
Total volume:	7.35 Litres
Weight when empty:	707g
Maximum recommended carry weight:	5 kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute

Product Information Sheet



- 1 x Laryngoscope Handle
- 1 x Blade Mac #4
- 1 x Blade Mac #3
- 1 x Blade Mac #2
- 1 x Syringe 10ml
- 1 x Artery Forceps
- 1 X Extension Tube
- 2 x Gloves
- 1 x Tape
- 3 x Lubricant
- 2 x Battery C size
- 1 x CO2 Detector
- 1 x Antimicrobial Filter
- 2 x Tube ties
- 1 x ODD (oesophageal detection device)



- 7 x ETT
- 1 x Magill forceps Adult
- 1 x Magill forceps Paed
- 2 x Stylet

Neann Three Quarter Airways Kit (Paediatric panel) 50-001-0642-0422



- 1 x Blade Mac #1
- 1 x Blade Mill #1
- 1 x Blade Mill #0
- 1 x Blade Mill #00
- 5 x ETT (Size 3, 3.5, 4, 4.5, 5)
- 2 x Stylet

Also has ETT sizing chart

