

Product Information Sheet



Product Code: 50-001-0652-0386

Product Name: NEANN Advanced Airway Kit

Purpose / Function / Design

Originally based of a design for NSW Ambulance Service to store and transport airway and intubation equipment

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Kit clearly named: Allows ease of identification
- Antimicrobial materials used: EnvironFrei improves infectious controls minimizing cross contamination
- aVector labelling : ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability, allows ease of cleaning
- Coated Webbing: Allows ease of cleaning and minimizes cross contamination
- Rigidity to lid and base : maintains form and adds impact resistance
- Heavy duty zip placed to give maximum access to interior

Internal


- MEDIS labelling: Improves on scene efficiencies through quick identification of contents
- Customized design: Eliminates “Search & Select” issues
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Paediatric ETT Size Chart included: Provides on scene efficiencies
- Internal flip panel folds out: all contents can be viewed at once
- Clear sleeve: Provides easy storage of documents

Specifications

External dimensions:	405mm (L) x 255mm (W) x 90mm (H)
Total volume:	9 Litres
Weight when empty:	967g
Maximum recommended carry weight:	3kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated Fabric
Internal structure:	Recycled Corflute

Date: Oct 2020

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365