

Product Information Sheet



Product Code: 50-001-0185-0637

Product Name: NEANN Transport Kit

Purpose / Function / Design

Easy to use, well-padded bag, designed for Intensive Care Units or Hospital in the Home to store and transport medical equipment and accessories.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions

Ease of cleaning, improves infection control minimizing cross contamination:

- Antimicrobial materials used: EnvironFrei
 - aVector labelling: ease of cleaning
 - Reinforced EnvironFrei fabric: Increased durability
 - Magnetised, replaceable and padded aVector handles
 - Coated webbing
- Padded, adjustable and replaceable shoulder sling: hands free carrying
 - Magnetic and replaceable sharps holder: secure holder for sharps waste, avoids contamination inside bag. Bag supplied with 600ml holder, option to purchase Larger Bio-can size
 - “Sharps attach” label: quickly re-attach Sharps holder, no search time **(Sharps Holder and Container are sold separately)**
 - Foam rigidity to outer walls, lid and base: maintains form and adds impact resistance
 - Heavy duty double zip with pull handle: placed to give maximum access to interior
 - Clear ID window: bag customization option with business card
 - Grab Handle on front: ease of moving the bag
 - External compartment with magnetic lid flap: auto close

Internal

- Water resistant anti-microbial fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Pocket (paper) labelling: Improves on scene efficiencies through quick identification of contents, can be personalized by end user
- Elasticated pockets on divider walls and wall of main compartment for stowing smaller items

Specifications

| | |
|-----------------------------------|--------------------------------------|
| External dimensions: | 440mm (L) x 280mm (W) x 250mm (H) |
| Total volume: | 30.8 Litres |
| Weight when empty: | 1.50kg |
| Maximum recommended carry weight: | 7kg |
| External fabric: | Medium weight reinforced EnvironFrei |
| Internal fabric: | Lightweight coated Fabric |
| Internal structure: | Foam |

Product Information Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Information Sheet



Shown with **Sharps container holder 600ml (2 mags) mag strap 50-001-3079-0386** & **Sharps Container 600ml (petal opening) 50-065-364151-0386** (separate products):



Shown with **NEANN Bio-can 1.4L Sharps container holder (2 mags) mag strap 50-001-8184-0386** (separate product):



Shown with **NEANN FITTANK 1.8L Sharps container holder (2 mags) mag strap 50-001-8192-0386** (separate product):



Shown with **NEANN 3/4 Airway 50-001-3079-0386** & **NEANN 3/4 Drug 50-065-364151-0386** (separate products):

