

Product Information Sheet



Product Code: 50-001-8733-0637

Product Name: OxyresQ Aquatic Kit & Pouch

Purpose / Function / Design

Designed to store and transport Oxygen cylinder, regulator, masks and related equipment in Aquatic conditions. Pouch will carry Oropharyngeal tubes, BVM and suction.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Anti-microbial materials used: EnvironFrei improves infectious controls minimizing cross contamination.
- aVector labelling: ease of cleaning and minimizes cross contamination
- Reinforced EnvironFrei fabric: Increased durability and allows ease of cleaning.
- Magnetised, replaceable and padded aVector handles: Minimizes cross contamination.
- customer logo included
- Foam rigidity to outer walls, lid and base : maintains form and adds impact resistance
- Heavy duty coated zip : prevents moisture ingress and ease of cleaning.
- Non-metal zip sliders : prevents corrosion in marine environments.
- Yellow zip pull tab extension: good visibility and ease of grab when wearing gloves
- Yellow grab handle at cylinder base ends.
- Adjustable backpack straps: keep hands free when carrying.
- Pouch has magnet attachment to Oxy boot for ease of attachment and movement.
- Keyhole from pouch to boot : to allow hose connection from cylinder to suction
- Clear windows both sides and top of pouch : to view gauge and pouch contents

Internal

- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- Water resistant fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Internal compartments have elasticated pockets inside for compartmentalised storage

Specifications

External dimensions:	590mm (L) x 150mm (W) x 320mm (H)
Total volume:	28 Litres
Weight when empty:	2.1 kg
External fabric:	Medium weight reinforced EnvironFrei
Internal fabric:	Lightweight coated fabric
Internal structure:	12mm Polyfoam

Product Information Sheet



Product Information Sheet



Product Information Sheet

