

## Compare Portable Devices by Features

|                          | Portable LED Devices                    |                             |                             |                          |  |  |
|--------------------------|---|-----------------------------|-----------------------------|--------------------------|--|--|
|                          | EMS PRO                                 | LED+                        | LEDX                        | PEDI2                    | NEO  | R  |
| Number of LEDs           | 26                                      | 28                          | 32                          | 14                       | 3  | 3  |
| Color of LEDs            | Orange (18)<br>+ Red (6) + White<br>(2) | Orange (22)<br>+ Red (6)    | Orange (24)<br>+ Red (8)    | Orange (10)<br>+ Red (4) | Orange (veins) +<br>Green (arteries)<br>+ White (exam) | Orange (veins) +<br>Green (arteries)<br>+ White (exam) |
| Color Selection          | Yes                                     | Yes                         | Yes                         | Yes                      | Yes  | Yes  |
| Brightness Selection     | No                                      | No                          | No                          | Yes                      | Yes, 3 Levels  | Yes, 3 Levels  |
| Ring Diameter (mm)       | 21                                      | 27                          | 31                          | 18                       | N/A  | N/A  |
| Battery                  | 2 x AA Alkaline                         | Rechargeable<br>Lithium Ion | Rechargeable<br>Lithium Ion | 2x AAA Alkaline          | 1 x AAA Alkaline                                       | 1 x AAA Alkaline                                       |
| Battery Run Time (mins)  | 550                                     | 350                         | 130 / 160<br>(32 / 24)      | -                        | -  | -  |
| Dimensions (mm)          | 116 x 55 x 23                           | 99 x 57 x 21                | 102 x 65 x 21               | 85 x 40 x 17             | 96.1 x 19 x 16.5                                       | 96.1 x 19 x 16.5                                       |
| Weight, inc. battery (g) | 108                                     | 77                          | 83                          | 64                       | 29   | 29   |
| Price                    | \$379                                   | \$549                       | \$679                       | \$249                    | \$299  | \$329  |