

Product Data Sheet



Product Code: 50-001-0656-0637

Product Name: NEANN BPR

Purpose / Function / Design

This Back Pack Resuscitator is designed to store and transport a Medical Oxygen cylinder (C, CC, CD (Inhalo) sizes) and related resuscitation equipment.

Features & Benefits

External

- Reflective tape: Provides high visibility in low light conditions
- Kit clearly labelled : Ease of identification

Ease of cleaning, improves infection control minimizing cross contamination:

- Antimicrobial materials used: EnvironFrei
- aVector labelling : ease of cleaning
- Reinforced EnvironFrei fabric: Increased durability

- Reinforced base: Abrasion resistant fabric
- Padded, low profile carry handles at top and side : user comfort, positioned so bag can be accessed in variable storage positions
- Padded, adjustable back pack straps with hip belt and cover : keep hands free when carrying, ergonomic weight distribution and easy stowage when not in use to reduce catch points
- Rigidity to outer walls, lid and base : maintains form and adds impact resistance
- Heavy duty zip placed to give maximum access to interior
- Clear ID window: bag customization option with business card

Internal

- MEDIS labelling: Improves on scene efficiencies through quick identification of contents
- Customized design: Eliminates "Search & Select" issues
- Everything has a place and is in its place: Provides a cost saving as eliminates overstocking
- Water resistant anti-microbial fabric: Allows ease of cleaning and minimizes cross contamination
- Lightly coloured fabric: Allows contamination to be easily identified
- Magnetic removable BVM pouches with locator labelling : dedicated location in bag
- Clear label sleeve on BVM Paediatric pouch and six pockets : labels can be added to customize for future proofing, improves on scene efficiencies through quick identification of contents, can be personalized by end user
- Clear lid on pouches and pockets : easily view contents
- Stretchy pockets on main compartment walls : stowing smaller items and suction consumables
- Shock absorber padding around Nevac variable suction device : protection
- Magnetic flaps on drug pouches : one hand usage and auto close
- Rotatable and removable suction canister holder
- Cylinder cradle with padded straps : security and protection against cylinder impact, holds both C sized and Inhalo cylinders

Specifications

External dimensions:

600mm (L) x 320mm (W) x 200mm (H)

Product Data Sheet



Total volume:	38 Litres
Weight when empty:	3.46 kg
Maximum recommended carry weight:	11 kg
External fabric:	Medium weight reinforced EnvironFrei
Staerksager:	Abrasion resistant fabric
Internal fabric:	Lightweight coated Fabric
	Foam and Recycled Corflute



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365

Product Data Sheet



Head Office
RAPP Australia Pty. Ltd.
160 McClelland Ave, PO Box 120,
Lara, Victoria, 3212, Australia
ABN 35 054 189 365